



## Laboratory Products

## Color System



Pressing has prepared the Color System Kit specifically in order to allow dental technicians to colour the surfaces of TSM Acetal Dental and acrylic resins according to their requirements.

The Colour System kit comprises:

- N° 16 types of pigments in pottery jars;
- N° 1 bottle of 10 ml of "Acelux"(methyl methacrylate based photosensitive lacquer)
- N° 1 mixing palette made of ceramics;
- N°1 bottle of 10 ml of "Lux Clea"(Degreasing agent and solvent for brushes);
- N°2 brushes (Type "00" and "1")

**Cod. 00300**

## Mini Color System



The Mini Color System kit comprise:

**Cod. 00299 V Type "Vita"**

N°10 types of pigments in pottery jars

**Cod. 00299 I Type "Ivoclar"**

N° 9 types of pigments in pottery jars

**Cod. 00297 Type "Special USA ( Doctor Kit )"**

N° 6 types of pigments in pottery jars

N°01 flacone da 10 ml "Acelux"

N° 01 flacone da 10 ml di Lux Clea

Pigments code:

<b>Cod. 00310</b> White	<b>Cod. 00311</b> Black
<b>Cod. 00314</b> Cyan Blue	<b>Cod. 00315</b> Purple
<b>Cod. 00318</b> Neutral brown	<b>Cod. 00319</b> Brown red
<b>Cod. 00322</b> Transluced grey blue	<b>Cod. 00323</b> Ochra
<b>Cod. 00326</b> "USA -A"	<b>Cod. 00327</b> "USA -B"
<b>Cod. 00312</b> Yellow	<b>Cod. 00313</b> Rosso Magenta
<b>Cod. 00316</b> Green	<b>Cod. 00317</b> Orange
<b>Cod. 00320</b> Opaque grey	<b>Cod. 00321</b> Transluced grey
<b>Cod. 00324</b> Orange red	<b>Cod. 00325</b> Purple red
<b>Cod. 00328</b> "USA -C"	<b>Cod. 00329</b> "USA -D"

## Fast Tray



Fast Tray is a photopolymerising material available in 50 pieces packages in 3 colours, blue, neutral and pink. Fast Tray is very rigid once cured and is used in the construction of special trays and occlusal registration bases. It has an even thickness which does not need to be mixed and is easily shaped, any excess material can be used again if stored correctly (e.g. in the box provided, away from UVA rays) Polymerization is achieved when using UVA type lamps.

**Cod. 14001** Neutral

**Cod. 14002** Pink

**Cod. 14003** Blue

## Gyflux



Photopolymerizable sealing film for models (applied on plaster)

**Cod. 00090** Contents 30 cc

**Cod. 00091** Contents 50 cc

## Acecril



Specific product for the chemical bonding of TSM Acetal Dental and acrylic resins

**Cod. 00082**

## Accessories



Oral tips, dispenser

## Pre-Mark



CFC - free Spray for occlusion

Content 75 ml

**Cod. 0205** Green

**Cod. 0206** Red

## Gingi-sil

Hardness 85 shore A  
Addition Cured Hydrophilic Vinyl  
Polysiloxane for the  
development of removable  
gingival masks.



Packaging 2 tube of 48 ml + 12  
mixer

**Cod. 12380**

## Acelux



Acelux is a photosensitive light cure resin.

It is used as a transparent light curing stain lacquer for sealing stains on TSM Acetal Dental and acrylic resin in conjunction with the Pressing Colour System. Acelux needs to be photopolymerised using UVA lamps in order to maximise its characteristics.

We recommend photopolymerised lamps are used, which produce a polymerisation spectrum from 320 to 550 Nm. (e.g. using the Pressing Complex Lux S-8 (light cure unit))

**Cod. 00302/1** content 2 x 10 ml

**Cod. 00095** content 30 ml

## Sil-Press



Sil Press is a two-component silicone elastomer which cures at room temperature by polyaddition reaction.

Available in the following versions:

- Type 10 shore A content kg.1+1 Cod. **00352-10**
- Type 10 shore A content kg.5+5 Cod. **00353-10**
- Type 21 Shore A content kg.1+1 Cod. **00352-21**
- Type 21 Shore A content kg.5+5 Cod. **00353-21**
- Type 28 Shore A content kg.1+1 Cod. **00352-28**
- Type 28 Shore A content kg.5+5 Cod. **00353-28**

## Sil-Labor



Sil-labor is a silicone elastomer which cures at room temperature by polycondensation reaction, specific for templates.

Contents:

- content kg.1,5 + 1 catalyst type (Classic) 85 Shore A Cod. **00360**
- type (Hard) 95 Shore A Cod. **00370**
- content kg. 5 + 3 catalyst type (Classic) 85 Shore A Cod. **00361**
- type (Hard) 95 Shore A Cod. **00371**

## Gemini



Silicon addition 1:1

Gemini is a two component silicone elastomer which cures at room temperature by polyaddition reaction, specific for templates.

- |                 |         |               |                      |
|-----------------|---------|---------------|----------------------|
| Type 82 Shore A | kg. 1,5 | A + Kg. 1,5 B | <b>Cod. 00380-82</b> |
| Pink            | kg. 5   | A + Kg. 5 B   | <b>Cod. 00381-82</b> |
| Type 90 Shore A | kg. 1,5 | A + Kg. 1,5 B | <b>Cod. 00380</b>    |
| Light Blue      | kg. 5   | A + Kg. 5 B   | <b>Cod. 00381</b>    |

## Pressing AV Abutment



The UV Abutment is a photopolymerising space resin which requires the use of UVA lamp.  
 The product is applied directly on the die preps and once cured it is resilient to steam and abrasive forces. Each application is of about 15 microns in thickness. When applying the resin, various thicknesses can be achieved, depending on the type of work carried out. UV Abutment has a thinner that can be used to a maximum of 20% in order to dilute the viscosity of the materials.  
 UV Abutment comes in a choice of four colours

**Cod. 0211** Blue - **Cod. 0212** Red  
**Cod. 0213** Yellow - **Cod. 0214** Green

## Poli Corflex



Polishing felts for Corflex Orthodontic

Contents: 11 pieces

**Cod. 00085**

Pressing Dental Srl

Via Onesto Scavino, 4  
 Falciano (Zona industriale Rovereta)  
 47891 - Rep. di San Marino

Tel. dall'Italia  
 0549-909948 - FAX 0549-909958  
 From other countries  
 \*\*378-909948 - FAX \*\*378-909958



## Corflex Brilliant



Aromatized polishing fluid for soft resins

Contents: 25 cc

**Cod. 00080**

## Separating Pressing



The "Separating Pressing" is a special insulating material for:

Acrylic teeth / plaster  
Acrylic resin / Plaster  
Plaster / Plaster

Contents: 500 ml

**Cod. S-500**

## Expo Model



A simple solution for your exhibition moulds  
Epoxy Resin compound

Packaging: 500 gr. compound A + 100 gr. B

**Cod. E-My** Yellow

**Cod. E-Mr** Red

**Cod. E-MG** Green

**Cod. E-MB** Blue

## J100 Tube



Containers for resin injection for the following Pressing machines: Actio, J-100, J-200 Evolution

**Cod. 00040** Contents 60 tubes + 60 caps

**Cod. 00041** Contents 75 tubes (no caps)

## Soft Tube



Soft Tubes for the following injection machines: MG.Newpress, Polyapress etc.

Contents 60 tubes + 60 caps

**Cod. 00042**

## Tube Lubricant



Special lubricant suitable for high temperatures for the following Pressing machines: Actio, J-100, J-200

Contents : 50 cc.

**Cod. 00050**

Lubricant only for the following machines: MG.Newpress, Polyapress etc

Contents : 50 cc

**Cod. 00051**

## Universal Polish



Universal polishing paste for acetal and acrylic resins.

Contents : 50 cc.

**Cod. 00060**

## Super Polish



La nuova pasta lucidante Universale per la resina termoplastica F.J. P.

Utilizzabile anche su tutte le resine a base di nylon.

Conf. da 50 cc

**Cod. 00060-S**

## Corflex Polish



Silicone polishing gel for positioners and gum shields.

Contents: 50 cc.

**Cod. 00060**

## Marble Stone Evolution



Special plaster recommended to develop models for TSM Acetal Dental processing.

Contents: Kg. 4 + P.L.

**Cod. 00084-E**

## Odilux



Denti in acrilico a tre strati della ditta UNIDESA (Spagna) in colorazione "IVOCLAR"

Carta delle forme composta da 20 forme incisivi superiori, 10 forme incisivi inferiori, 7 forme molari superiori e inferiori.

## Odipal



Denti in acrilico a tre strati della ditta UNIDESA (Spagna) in colorazione "VITA"

Carta delle forme composta da 20 forme incisivi superiori, 10 forme incisivi inferiori, 12 forme molari superiori e inferiori.

### Odipal PR

Versione solo posteriore con un profilo vestibolare più allungato e caratterizzato, senza un profilo occlusale definito, come nei denti naturali. Carta delle forme composta da 3 forme.

## Masky



The mask for your metal; it is a thermosetting product created to mask the metallic parts of prostheses.

Available:

**Code M-0100** White 6-8 gr.( n.b. is possible to obtain other products mixing it with the colour system)

**Code M-0105R** Pink 3 gr.

**Code M-I (Masky Investment)** Coating for Masky Kg. 1.

## Fil wax



White Wax suitable for making channels during the processing.

Available diameters:

**Cod. 02160** 2 mm gr. 140

**Cod. 02161** 3 mm gr.160

**Cod. 02163** 5 mm gr.160

**Cod. 02162** 4 mm gr. 160

**Cod. 02164** 9 mm gr.160