



Corflex Orthodontic

Classificazione secondo la direttiva 93/42CEE
Classe II A CE 0546

Corflex-Orthodontic is a synthesis product obtained from a mixture of high molecular weight copolymers of ethylene and vinyl acetate. The association of its components results from the research efforts carried out by Pressing laboratories; these enabled Pressing to produce a thermoplastic material with the ideal features to meet a whole range of requirements, so much so that Corflex-orthodontic has become the first choice product of its kind. Rapid and recurrent mechanic stresses do not cause any cracking in the Corflex-orthodontic; it displays a high degree of flexibility and a rubber-like performance. Its very low water absorption and excellent resistance to acids are further noteworthy features.

The Corflex-orthodontic is available in a wide colour range, it is so transparent that controls can be carried out even after its placement in the oral cavity. Applications: Corflex-orthodontic is widely used to make orthodontic positioning devices, oral devices and gum shields for sport activities.

Corflex Orthodontic underwent the following tests:

Citotoxicity test (EN 30093-part 5)

Cell vitality assessment by means of propidium iodide

Cell vitality assessment by means of neutral red absorption

Genotoxicity test (EN 30093-part 5)

Mutagenicity assessment on Ames Typhimurium Salmonella

Chromosomic aberration assessment

Assessment of exchanges between chromatid brothers

Test di Anallergicità (ISO/DIS 10093-part 10)

Cutaneous irritation tests

Maximum sensitisation test

Water absorption and solubility in water (ISO 10477)

Products codes:

Cod. 00997 Colour shade

Cod. 01000 Assorted colours gr. 600 (except white and black)

Cod. 01001 Transparent (Packaging of gr. 400)

Cod. 01002 White **Cod. 01003** Yellow **Cod. 01004** Green

Cod. 01005 Pink **Cod. 01006** Red **Cod. 01007** Fuchsia

Cod. 01008 Light Blue **Cod. 01009** Blue **Cod. 01010** Black



Gum shields



Poly Corflex



Polishing felts for Corflex Orthodontic

Contents: 11 pieces

Cod. 00085A Polishing felts 1st pass

Cod. 10085B Polishing felts 2nd pass

Corflex Brilliant



Aromatized polishing fluid for soft resins

Contents: 25 cc

Cod. 00080

Pressing Dental Srl

Via Onesto Scavino, 4
Falciano (Zona industriale Rovereta)
47891 - Rep. di San Marino

Tel. dall'Italia
0549-909948 - FAX 0549-909958
From other countries
**378-909948 - FAX **378-909958